



PSL

Portable Slit Lamp

A full featured slit lamp, in the palm of your hand.

PSL Portable Slit Lamp

The Reichert® **PSL, Portable Slit Lamp**, utilizes a high luminance, natural daylight LED, coupled with a lightweight, compact design, to give you the freedom of using a slit lamp in the palm of your hand. High quality optics provide a clear and detailed view of the eye.

Small and Lightweight

The PSL Slit Lamp weighs about 1.5 lbs. (0.68 kg) including battery and fits easily and comfortably in your hand. The ergonomic, portable design and convenient carrying case make it an ideal instrument for non-ambulatory patients, nursing home visits, and veterinary work.

High Luminance Natural Daylight LED

The PSL's high luminance light source provides about 10,000 lux of natural daylight illumination and has a lifespan of 50,000 hours or more; providing bright, consistent light, as well as the convenience of not having to replace bulbs.

Long Battery Life

The PSL Portable Slit Lamp comes with two long-lasting batteries that can power the unit at maximum illumination for up to two hours. For added convenience, the PSL also comes with its own battery charger.

Specifications

Catalog Number	15090 - PSL Portable Slit Lamp
Magnification	10x & 16x eyepieces included
Diopter Adjustment	± 7D
Visual Field	Dia. 0.39 in. (10 mm)
PD Adjustment Range	10x eyepiece: 1.97 - 2.76 in. (50 - 70 mm) 16x eyepiece: 2.17 - 2.76 in. (55 - 70 mm)
Working Distance	3.15 in. (80 mm)
Slit Width	0 - 0.43 in. (0 - 11 mm)
Slit Length	0.43 in. (11 mm)
Filters	Cobalt Blue, Red Free, Color Temperature Conversion
Slit Rotation	± 30°
Light Source	High Luminance Natural Daylight LED
Battery	Li-ion battery 7.4V, 680 mAh (2 hrs. operation at max. illumination)
Charger	Input: 100-240V AC, 50-60 Hz Output: 12V DC, 500 mA
Standard Accessories	Carrying case, Battery (2), Battery charger (1), Forehead support (1), Diopter adjustment bar (1)
Dimensions (Slit Lamp)	4 W x 8.38 H x 6.25 D in. (102 x 213 x 159 mm)
Weight	Approx. 1.5 lbs. (680 g)



Visit www.reichert.com or contact your Authorized Reichert Technologies Distributor.

Important Notice: For your own protection, purchase instruments from and have them serviced only by an Reichert Technologies or an Authorized Reichert Technologies Distributor.

15090-110-Rev-A ©2012 AMETEK, Inc. All rights reserved. - Specifications subject to change without notice. Reichert is a registered trademark of Reichert, Inc. Reichert Technologies and "Advancing Eye Care. Preserving Sight." are trademarks of Reichert, Inc. AMETEK is a registered trademark of AMETEK, Inc.



Reichert, Inc.
3362 Walden Avenue, Depew, NY 14043
Toll Free: 888-849-8955 • Tel: (716) 686-4500
Fax: (716) 686-4545 • www.reichert.com

ISO 9001 Certified
Reichert products are designed and manufactured under quality processes meeting ISO 9001 requirements.