

Specifications

Microscope

Type	Stereoscopic Microscope		
Magnification selection	2 steps by Objective Lens rotation		
Objective	1X, 1.6X		
Eyepiece	10X, 16X (16X Eyepiece is optionally available in some regions.)		
Magnification ratio	Objective	Eyepiece	Magnification (Field of View)
		1X	10X (18mm dia.)
	1.6X	16X	16X (14.5mm dia.)
		10X	16X (11.25mm dia.)
Pupillary adjustment	16X	25.6X (9mm dia.)	
	10X Eyepiece	55mm-82mm	
Diopter adjustment	16X Eyepiece	51mm-78mm	
	10X Eyepiece	±8D	
	16X Eyepiece	±10D	

Slit Illumination

Slit width	Continuous from 0 to 9mm (Circle at 9mm)		
Slit Length	Continuous from 1 to 8mm		
Aperture diameter	9, 8, 5, 3, 2, 1, 0.2mm		
Slit angle	0-180°		
Slit inclination	5°, 10°, 15°, 20°		
Filters	Cobalt blue, Red-free, 13% ND, Heat absorbing, UV Cut (Normal use), IR Cut (Normal use)		
Illumination lamp	6V 20W Halogen		

Base

Longitudinal movement	90mm
Lateral movement	100mm
Fine base movement	15mm
Vertical movement	30mm

Chinrest

Vertical movement	80mm
Fixation target	6V 0.2A Tungsten lamp

Power

Input (Primary)	AC 120V/220V/230V/240V
Frequency	50-60Hz
Power Consumption	30VA

Dimension & Weight

Table type	550mm(W) x 430mm(D) x 770-800mm(H), 21Kg
Unit type	440mm(W) x 410mm(D) x 770-800mm(H), 20Kg

SL-3G

Slit Lamp



* Not available in all countries, please check with your distributor for availability in your country
 * Subject to change in design and/or specifications without advanced notice.

IMPORTANT In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.

TOPCON CORPORATION TOPCON CORPORATION TOPCON CORPORATION

TOPCON CORPORATION

TOPCON MEDICAL SYSTEMS, INC.

111 Bauer Drive, Oakland, NJ 07436, U.S.A.
 Phone: +1-201-599-5100 Fax: +1-201-599-5250 www.topconmedical.com

TOPCON CANADA INC.

110 Provencher Avenue, Boisbriand, QC J7G 1N1 CANADA
 Phone: +1-450-430-7771 Fax: +1-450-430-6457 www.topcon.ca

TOPCON EUROPE MEDICAL B.V.

Essebaan 11, 2908 LJ Capelle a/d IJssel, P.O.Box 45,
 2900 AC Capelle a/d IJssel, THE NETHERLANDS
 Phone: +31-(0)10-4585077 Fax: +31-(0)10-4585045
 E-mail: medical@topcon.eu, www.topcon-medical.eu

Italy Office: Viale dell'Industria 60, 20037 Paderno Dugnano, (Milano), ITALY
 Phone: +39-02-9186671 Fax: +39-02-91081091 E-mail: info@topcon.it, www.topcon.it

DANMARK OFFICE

Praestemarksvej 25, 4000 Roskilde, DANMARK
 Phone: +45-46-327500 Fax: +45-46-327555
 E-mail: topcon@topcondanmark.dk www.topcondanmark.dk

IRELAND OFFICE

Unit 276, Blanchardstown, Corporate Park 2 Ballycoolin Dublin 15, IRELAND
 Phone: +353-18975900 Fax: +353-18293915 E-mail: medical@topcon.ie, www.topcon.ie

TOPCON S.A.R.L.

BAT A1 3 route de la Révolette 93206 SAINT DENIS CEDEX
 Phone: +33-1-49-21-23-23 Fax: +33-1-49-21-23-24 E-mail: topcon@topcon.fr, www.topcon.fr

TOPCON DEUTSCHLAND MEDICAL G.m.b.H.

Hanns-Martin-Schleyer Strasse 41, D-47877 Willich, GERMANY
 Phone: +49-(0)2154-8850 Fax: +49-(0)2154-885177 E-mail: med@topcon.de, www.topcon.de

TOPCON SCANDINAVIA A.B.

Neogatan 2, P.O.Box 25, 43151 Mölndal, SWEDEN
 Phone: +46-(0)31-7109200 Fax: +46-(0)31-7109249 E-mail: medical@topcon.se, www.topcon.se

TOPCON ESPAÑA S.A.

HEAD OFFICE
 Frederic Mompou 4 Esc. A Bajos 3, 08960 Sant Just Desvern Barcelona, SPAIN
 Phone: +34-93-4734057 Fax: +34-93-4733932 E-mail: medica@topcon.es, www.topcon.es

TOPCON (GREAT BRITAIN) LTD.

Topcon House, Kennet Side, Bone Lane, Newbury, Berkshire RG14 5PX United Kingdom
 Phone: +44-(0)1635-551120 Fax: +44-(0)1635-551170 E-mail: medical@topcon.co.uk, www.topcon.co.uk

TOPCON POLSKA Sp. z o. o.

ul. Warszawska 23, 42-470 Siewierz, POLAND
 Phone: +48-(0)32-6705045 Fax: +48-(0)32-6713405 www.topcon-polska.pl

TOPCON SINGAPORE MEDICAL PTE. LTD.

1, JALAN KILANG TIMOR #09-01 PACIFIC TECH CENTRE SINGAPORE 159303
 Phone: +65-68720606 Fax: +65-67736150
 E-mail: medical_sales@topcon.com.sg www.topcon.com.sg

REPRESENTATIVE OFFICE IN INDONESIA

Level 38, Tower A, Kota Kasablanka unit GH-04, Jl. Casablanca, Kav 88, Jakarta, Indonesia 12870
 Phone: +62-21-2963-8004

TOPCON INSTRUMENTS (MALAYSIA) SDN.BHD.

No. D1, (Ground Floor), Jalan Excella 2, Off Jalan Ampang Putra,
 Taman Ampang Hilir, 55300 Kuala Lumpur, MALAYSIA
 Phone: +60-(0)3-42709866 Fax: +60-(0)3-42709766

TOPCON INSTRUMENTS (THAILAND) CO.,LTD.

77/62 Sinsathorn Tower, 37th Floor, Krungthornburi Rd., Klongtonnai,
 Klongsarn, Bangkok 10600, THAILAND
 Phone: +66-(0)2-440-1152-7 Fax: +66-(0)2-440-1158

TOPCON CORPORATION BEIJING OFFICE

Block No.9, Kangding Street Beijing Economic-Technological Development Area,
 Beijing, 100176, CHINA
 Phone: +86-(0)10-6780-2799 Fax: +86-(0)10-6780-2790

PERFORMANCE
 YOU CAN COUNT ON



A Classic Slit Lamp with Unparalleled Optical Quality

The new SL-3G, with a style that became the benchmark of slit lamps worldwide, combines classic versatility with the smoothest mechanism in its type.



Slit Lamp
SL-3G



Optical System

Topcon SL-3G's stereoscopic converging microscope with the use of innovative optics provides the practitioner with unsurpassed image clarity and sharpness. This reduces fatigue while providing a clearer field of view. Magnifications can easily and smoothly be changed between 1X to 1.6X using the conveniently located lever. Final magnification ratios are equivalent to 10X and 16X with the 10X eyepiece; 16X and 25.6X with the 16X eyepiece.

Bright, Even, Consistent Illumination

With the use of a powerful halogen bulb the SL-3G produces a consistent, steady lighting for observation. Its 6V 20W halogen bulb resembles natural daylight illumination for optimum color rendition. The intense beam allows for narrower slit apertures and finer detail observation. Longer lamp duration and a steady brightness during the lamp's life are additional advantages.

Smooth, Effortless Controls

The SL-3G features the Topcon Omni-Directional control lever that provides the smoothest movement for precise, fatigue-free adjustments.



Slit Beam Inclination

Added versatility in gonioscopy and fundus observation is provided by the capability to incline the illumination unit up to 20° from the observation axis. Horizontal sections can also be observed with this feature.

Four Built-in Filters

The SL-3G comes with four built-in filters into the illumination column to facilitate special examinations. The filters are easily inserted into the illumination path by a quick change lever. The filters include: Cobalt blue, Red-free, 13% Neutral Density and Heat Absorbing.

Continuous Adjustment of Slit Aperture

Both slit width and length are continuously adjustable. The width can be adjusted from 0 to 9mm and the length from 1 to 8mm. A rotating collar allows for quick adjustments of 0.2, 1, 2, 3, 5, 8 and 9mm.

Precise Angular Measurements

An angular measurement scale is engraved on the SL-3G housing to simplify the fitting of astigmatic contact lenses.



Optional Accessories

- Hruby lens
- 10X measuring eye piece
- The SL-3G accepts applanation tonometers models R900 and T900
- The SL-3G can be integrated with the Topcon Laser Photocoagulator LC-300G by using the SA-1G Slit Lamp Delivery System

